

SALINE MONITORING SYSTEM.....MAKING MONITORING EASY



1

Attach the Saline bottle apparatus to the setup of IR sensors, place the drip chamber at correct positions across the IR sensors properly.

2

After the setup is attached the system will start working after giving power supply.

3


The first IR sensor will start detecting the drop and estimate the time to get the whole saline bottle empty.

4

The first alert message will be sent to the user, after the drip chamber getting empty a second alert will be sent as a precautionary alert.

5

After getting the alerts, the concerned medical authority can come and change or remove the saline bottle and reset the system.



**SALINE MONITORING SYSTEM HELPS
AND EASES THE WORK OF DOCTORS
AND NURSES GIVING ACCURATE
RESULTS**

